

2X YourTaq™ PCR Master Mix

LOT: See product label EXPIRY DATE: See product label

ORDERING INFORMATION

CAT. NO.	SIZE	PACKAGE CONTENT
BR0102201	200 reactions of 50 μl	4 × 1.25 ml YourTaq PCR Master Mix
BR0102202	1000 reactions of 50 µl	20 × 1.25 ml YourTaq PCR Master Mix

COMPONENT	COMPOSITION
YourTaq PCR Master Mix	Optimized 2X YourTaq PCR Master Mix
STORAGE	-20°C (until expiry date – see product label)

FEATURES

- Exceptionally pure YourTag DNA Polymerase
- · Optimized Master Mix for increased yield of amplification
- · Resistant to PCR inhibitor carry-over
- Excellent PCR specificity and sensitivity for a broad range of amplicons

APPLICATIONS

- · Routine and demanding PCR amplification up to 3 kb
- Suitable for amplification of low target copy number
- TA cloning

2X YourTag™ PCR Master Mix

DESCRIPTION

biotechrabbit™ YourTaq PCR Master Mix is optimized for high yield of amplification of 0.1–3 kb DNA targets, even from low copy number. The mix shows excellent PCR specificity and sensitivity for a broad range of amplicons. YourTaq PCR Master Mix is resistant to PCR inhibitors, such as blood (up to 20%), Ethanol or humic acid enabling PCR amplification from DNA templates with carry-over of PCR-inhibitors.

The 2X YourTaq PCR Master Mix contains pure biotechrabbit YourTaq DNA Polymerase, extremely high-quality dNTPs and optimized PCR buffer; thus, only template, PCR primers and PCR-grade water are added. In addition, the mix is suitable for amplification of GC-rich templates (up to 70%) pairing with 5X PCR Enhancer (BR1900201).

Info: Recommended annealing temperature is 2°C above primer Tm (use gradient PCR to check).

PROTOCOL

Prevention of PCR contamination

When assembling the amplification reactions, care should be taken to eliminate the possibility of contamination with undesired DNA.

- Use separate clean areas for preparation of samples and reaction mixtures and for cycling.
- Wear fresh gloves. Use sterile tubes and pipette tips with aerosol filters for PCR setup.
- Use only water and reagents that are free of DNA and nucleases.
- With every PCR setup, perform a contamination control reaction that does not include template DNA.

Standard PCR setup

The standard PCR protocol using biotechrabbit reaction buffer provides excellent results for most applications. Optimization might be necessary for certain conditions, such as the amplification of long targets, high GC or AT content, strong template secondary structures or insufficient template purity. In such cases, optimization of template purification (see biotechrabbit nucleic acid purification kits), primer design and annealing temperature is recommended.

The best conditions for each primer-template can be optimized with the following:

- Choosing the optimal quantities of template and primers
- Optimizing cycling conditions

BASIC PROTOCOL

- The Master Mix is designed to be used without any optimization as it has all necessary reaction components in optimal amounts for successful PCR.
- 5X PCR Enhancer (BR1900201) can be used to apply the mix for GC-rich target amplification (up to 70% GC).
- Thaw on ice and mix all reagents well.
- · Keep all reagents and reactions on ice.
- Pipet the master mix into thin-walled 0.2 ml PCR tubes.
- Add template and primers separately if they are not used in all reactions.

COMPONENT	VOLUME	FINAL CONCENTRATION				
2X YourTaq PCR Master Mix	25 μΙ	1×				
5X PCR Enhancer (optional, see BR1900201)	10 μΙ	1×				
Forward primer	Variable	0.2–1 μM				
Reverse primer	Variable	0.2–1 μM				
Template DNA	Variable	10 pg–1 μg				
Use 0.01-1 ng for plasmid or phage DNA and 0.05-1 μg for genomic DNA						
Nuclease free water	Variable					
Total volume	50 μΙ					

- For total reaction volumes other than 50 µl, scale reagents proportionally.
- Mix and centrifuge briefly to collect the liquid in the bottom of the tube.
- Place in the PCR cycler.

CYCLING PROGRAM

STEP	TEMPERATURE	TIME	CYCLES		
Initial activation	95°C	2 min	1		
Denaturation	95°C	30 s	25–35		
Annealing*	(55-68°C)	15–30 s	25–35		
	*Recommended annealing temperature is 2°Cabove Tm of primers.				
	Use gradient PCR to optimize the annealing temperature.				
Extension	72°C	30-60 s/kb	25–35		
Final extension	72°C	5 min	1		
	To extend all incomplete PCR products				
Storage in the cycler	4°C	Indefinitely	1		

• Add loading dye solution (see 6X DNA Loading Dye, BR0800301) to the reactions to analyze PCR products on a gel or store them at -20°C.

CERTIFICATE OF ANALYSIS

Quality Control

Functional assav

Human genomic DNA was amplified using the YourTag PCR Master Mix and specific primers to produce a distinct band.

Quality confirmed by: Head of Quality Control

SAFFTY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter: http://www.biotechrabbit.com/support/documentation.html.

USEFUL HINTS

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

CONTACT BIOTECHRABBIT

biotechrabbit GmbH

Volmerstr. 9a info@biotechrabbit.com 12489 Berlin. support@biotechrabbit.com Phone: +49 30 555 7821-10 Germany

www.biotechrabbit.com Fax: +49 30 555 7821-99



Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.

Trademarks: biotechrabbit™, YourTaq™ (biotechrabbit GmbH). valid from 05.02.2021

