

## cDNA Synthesis Kit

LOT: See product label

EXPIRY DATE: See product label

## ORDERING INFORMATION

CAT. NO.	SIZE	PACKAGE CONTENT
BR0400401	10 rxn of 20 $\mu$ l	10 $\mu$ l RevertUP™ II Reverse Transcriptase 100 $\mu$ l 5 $\times$ Reverse Transcriptase Buffer 25 $\mu$ l dNTP Mix (10 mM each) 5 $\mu$ l RNase Inhibitor 10 $\mu$ l Hexamer Primer 5 $\mu$ l Oligo (dT) Primer 1.5 ml PCR Grade Water
BR0400403	125 rxn of 20 $\mu$ l	125 $\mu$ l RevertUP™ II Reverse Transcriptase 1 ml 5 $\times$ Reverse Transcriptase Buffer 250 $\mu$ l dNTP Mix (10 mM each) 62.5 $\mu$ l RNase Inhibitor 125 $\mu$ l Hexamer Primer 62.5 $\mu$ l Oligo (dT) Primer 5 $\times$ 1.5 ml PCR Grade Water

COMPONENT	COMPOSITION
RevertUP II Reverse Transcriptase	RevertUP II Reverse Transcriptase, 200 U/ $\mu$ l in Storage buffer, containing 50% glycerol
5 $\times$ Reverse Transcriptase Buffer	Optimized 5 $\times$ Reverse Transcriptase buffer for cDNA synthesis
dNTP Mix (10 mM each)	Aqueous solution (pH 7.0) containing 10 mM each: dATP, dCTP, dGTP, dTTP sodium salts
RNase Inhibitor	RNase Inhibitor, 40 U/ $\mu$ l, in Storage buffer, containing 50% glycerol
Hexamer Primer	25 $\mu$ M Random Hexamer Primer
Oligo (dT) Primer	10 $\mu$ M Oligo (dT) Primer
PCR Grade Water	Ultrapure, sterile filtrated water, DNase-, RNase- and protease-free

STORAGE -20°C (until expiry date – see product label)

# cDNA Synthesis Kit

---

## FEATURES

- Highly efficient synthesis of long cDNAs ( $\geq 19$  kb)
- Excellent yields at temperatures up to 55°C
- High sensitivity reverse transcription from low abundance template
- Superior performance in demanding applications, including templates with a high degree of secondary structure

## APPLICATIONS

- First-strand cDNA synthesis for RT-PCR and qPCR
- Gene expression profiling
- RNA labeling and primer extension
- cDNA library construction

## DESCRIPTION

biotechrabbit™ cDNA Synthesis Kit provides superior components that ensure efficient first-strand cDNA synthesis from mRNA or total RNA templates. biotechrabbit RevertUP™ II Reverse Transcriptase enables highly efficient reverse transcription with increased thermostability. biotechrabbit RNase Inhibitor is a potent non-competitive inhibitor of RNases. The combination of highly efficient cDNA synthesis, effective RNase inhibition and pure dNTPs allows high yields of cDNAs of more than 19 kb.

For greater application flexibility, hexamer primers, allowing all RNAs in the reaction to be used as templates, and an oligo (dT) primer, for the synthesis of cDNA from only poly(A) tailed mRNA, are included.

## PROTOCOL

### Prevention of cDNA synthesis reaction contamination

RNase contamination is an exceptional concern when working with RNA. RNase A, providing most threat to RNA integrity, is a highly stable contaminant of any laboratory. To prevent RNA from degradation and to minimize possibility of contamination during cDNA synthesis; follow the guidelines below:

- Use separate clean areas for preparation of the samples and the reaction mixture.
- DEPC-treat all tubes and pipette tips or use certified nuclease-free labware with aerosol filters.
- Wear fresh gloves when handling RNA and all reagents.
- Always assess the integrity of RNA prior to cDNA synthesis in denaturing agarose gel electrophoresis.
- Use RNase free water and other reagents.
- To prevent RNA from degradation, add Ribonuclease inhibitor (optional) in to the cDNA synthesis reaction (20 units for 20  $\mu$ l reaction).

## Typical cDNA synthesis reaction set up

- Thaw on ice and mix very well all reagents.
- Assemble and keep all reactions on ice.
- To use time and reagents effectively, always prepare master mix for multiple reactions. For a master mix volume, always calculate the number reactions that you need plus one additional.
- Combine the following in an RNase-free reaction tube:

COMPONENT	VOLUME	FINAL CONCENTRATION
dNTP Mix (10 mM each)	2 $\mu$ l	1 mM (each dNTP)
RNase Inhibitor, 40 U/ $\mu$ l (optional)	0.5 $\mu$ l	1 U/ $\mu$ l
<i>Oligo (dT)<sub>12-18</sub> (10 <math>\mu</math>M) - or</i>	0.5 $\mu$ l	0.25 $\mu$ M
<i>Hexamer Primer (25 <math>\mu</math>M) - or</i>	1 $\mu$ l	1.25 $\mu$ M
<i>Gene Specific Primer (10 <math>\mu</math>M)</i>	0.5 $\mu$ l	0.25 $\mu$ M
5 $\times$ Reverse Transcriptase Buffer	4 $\mu$ l	1 $\times$
RNA Template	0.1–1 $\mu$ g total RNA or 50–500 ng mRNA (polyA)	
RevertUP™ II Reverse Transcriptase	1 $\mu$ l	10 U/ $\mu$ l
PCR Grade Water	Variable	
Total volume	20 $\mu$ l	

- Mix and collect the drops by centrifuging briefly.
- When using
  - Hexamer Primer, incubate 10 minutes at 30°C followed by 50–55°C for 20–60 minutes
  - Oligo (dT) or gene-specific Primer incubate at 50–55°C for 20–60 minutes.
- Inactivate enzyme at 99°C for 5 minutes.
- Collect the drops by spinning briefly.
- Store products at –20°C or proceed to next step, like PCR or qPCR.
- Use maximum 10  $\mu$ l of the cDNA synthesis reaction mix for PCR in 50  $\mu$ l volume.

# cDNA Synthesis Kit

---

## CERTIFICATE OF ANALYSIS

Quality Control

Functional Assay:

cDNA synthesis with specific primers, followed by quantitative PCR.

Quality confirmed by: Head of Quality Control

## SAFETY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter:

<http://www.biotechrabbit.com/support/documentation.html>.

## USEFUL HINTS

- Visit Applications at [www.biotechrabbit.com](http://www.biotechrabbit.com) for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

## CONTACT BIOTECHRABBIT

biotechrabbit GmbH  
Volmerstr. 9a  
12489 Berlin,  
Germany

[info@biotechrabbit.com](mailto:info@biotechrabbit.com)  
[support@biotechrabbit.com](mailto:support@biotechrabbit.com)  
[www.biotechrabbit.com](http://www.biotechrabbit.com)

Phone: +49 30 555 7821-10  
Fax: +49 30 555 7821-99



### *Legal Disclaimer and Product Use Limitation*

*Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.*

Trademarks: *biotechrabbit™, RevertUP™ (biotechrabbit GmbH).*

*valid from 25.04.2019*