




5,000 µl Biohit-Style-XL Tip (non-sterile)

SKU E1009-7001



Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Sterile	No
Design Tips	Biohit-Style-XL
Graduations	No
Volume	5000 µl
Pack Size	1 Bag x 250 Tips
For use with	Abimed: Discovery Comfort, LABMATE Optima Biohit®: eLine, ePet, mLine, Picus, Proline Eppendorf®: Reference 2, Research plus, Research pro, Xplorer Gilson®: Pipetman Concept, Pipetman P, Pipetman Ultra Hamilton®: SofTouch HTL: Discovery Comfort, Labmate Pro, Optipette Rainin®: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus Sartorius®: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta Thermo Fisher Scientific Inc.: Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus,

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 5,000 µl Biohit-Style-XL Tip (non-sterile)

Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de