



50 ml Centrifuge Tube, Skirted, Loose (Sterile)

SKU E1450-0400



Product Highlights

- Clear, sterile polypropylene tubes with green HDPE screw caps
- White writing surface and clear, printed graduations
- Perfect for centrifugation of cells: USP Class VI, non-cytotoxic and sterile (SAL 10^{-6}) material
- Certified RNase, DNase, DNA, heavy-metal and Pyrogen free
- Max. centrifugation rate: 6,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Tube Base	Skirted
Dimensions (W × H)	34.9 mm × 115.5 mm
g-Force max.	6000
Color	Clear
Volume	50 ml
Volume range	50 ml
Pack Size	500 Tubes (20 Bags × 25 Tubes)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

STARLAB International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 50 ml Centrifuge Tube, Skirted, Loose (Sterile)

STARLAB's new centrifuge tubes are designed for the centrifugation of cells, but can also be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. The sterility of assurance level (SAL) of 10^{-6} fulfills highest demands on sterility, meeting EN 556-1 standard (Sterilisation of medical devices). Racked tubes are delivered in recyclable cardboard racks.

Maximum centrifugation rates: Conical Tubes 15 & 50 ml: 12,000 x g and Skirted Tube 50 ml: 6,000 x g