

200 μ l UltraPoint[®] Yellow Graduated Tip, Loose (non-sterile)

TipOne[®] is the original refillable, reusable and recyclable pipette tip system. Recent improvements to the design of the system, based on customer-feedback, is why TipOne[®] is one of the leading tips for research labs.

Convenience and quality are at the heart of the system and TipOne[®] remains the same tip that you have used and trusted over the years.

Tip fit

One tip fits all! The tip collar is engineered to fit all popular pipettes.

Manufacture

The use of multiple cavity tools for the injection moulding process guarantees high conformity within one batch. The automated moulding and assembling procedures are carried out under clean-room conditions to eliminate contamination.

Unique collar

Our newest tips can be recognised by the jagged design on the tip collar.

Material

All TipOne[®] tips are manufactured from pure, virgin polypropylene to ensure the highest quality product.

Green tip design

A reduction of up to 40 % in the material used for the production of TipOne[®] tips helps to save valuable resources and is beneficial for the environment.

Highly polished moulds

For an extremely smooth surface to reduce sample retention.

Quality Control

Automated Quality Control procedures ensure that TipOne[®] is inspected for 100% performance. Every TipOne[®] product is certified RNase, DNase, DNA and Pyrogen free.

Loose tips are supplied in a convenient, self-standing bag which provides easy access and storage. The tamper-proof perforation strip on the bag enables easy opening. A resealable zip-lock provides contamination-free storage of opened product. The bag is made from recyclable polyethylene.


UltraPoint orifice

Thin-wall orifice helps prevent sample hang up and makes it easier to touch off the last drop of sample



Product Features

- **Graduated** at 10 µl, 50 µl and 100 µl
- **UltraPoint** thin-wall orifice helps prevent sample hang up and makes it easier to touch off the last drop of sample
- **Universal fit:** One tip fits all! The tip collar is engineered to fit all popular pipettes.
- **Highly polished moulds:** For an extremely smooth surface to reduce sample retention.
- **Material:** All TipOne® tips are manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified:** Every TipOne® product is certified RNase, DNase, DNA and Pyrogen free.

200 µl UltraPoint® Yellow Graduated Tip, Loose (non-sterile)			
	Description	Pack Size	Cat. No.
	200 µl UltraPoint® Yellow Graduated Tip, Loose (non-sterile)	1 Bag × 1000 Tips	S1113-1006

Technical Specifications	
Autoclavable	Yes
Color Code Wafer	Yellow
Material	Polypropylene
Pack Size	1 Bag × 1000 Tips
Sterile	No
Color	Yellow
Volume	200 µl
For use with	Abimed®, Biohit®, BRAND®, Capp®, Eppendorf®, Gilson®, Hamilton®, HTL, Nichiryo®, Rainin®, Sartorius®, Socorex®, STARLAB®, Thermo Scientific™